



# QUICK START GUIDE

**IMPORTANT!** Always read the operating instructions before starting up the unit to obtain full safety and operating instructions.

- 1** Place the unit on a solid and flat surface in a levelled position. Make sure that the room is ventilated and dry. Ensure air circulation around the unit and protect it against direct sunlight.
- 2** The door keys are affixed to the power cable on the back side of the unit.
- 3** Remove the protective film from the display.
- 4** Make sure that the plug is mounted correctly, the power supply is stable and that the unit is grounded and plugged to the correct voltage in accordance with the serial plate.
- 5** Press the on/off button to start the unit. We recommend letting the unit run for 24 hours before usage. Until the set temperature is reached, different alarms may occur.
- 6** A low battery alarm  may occur after switching on the unit. The battery will automatically recharge. Once the battery is fully recharged, the alarm icon will automatically disappear.
- 7** Do not load any items into the unit until the set temperature is reached.
- 8** Do not change the factory settings of the automatic defrost cycle as this can impair the performance of the unit (does not apply for -86°C ULT freezers).
- 9** External probes must be mounted through a porthole. Do NOT mount external probes or sensors through the door opening. This will affect the performance and stability of the unit and may cause condensation and/or ice build up on the door gasket and inside the unit.

To accept an alarm:



Press ENTER, then press the RIGHT KEY to highlight the OK key and then press ENTER once more to accept the alarm. Press ESC to return to the main screen.

Distributed By:

 **Azer Scientific**  
Serving Scientists Who Change The World\*

701 Hemlock Road,  
Morgantown PA 19543  
Tel: 1-877-770-2937  
Email: [Info@AzerSci.com](mailto:Info@AzerSci.com)  
[www.AzerSci.com](http://www.AzerSci.com)

5083030-05

**ARCTIKO**   
true cooling specialists