

## Tacta Mechanical Pipettes Perfectly Balanced



# Tacta Mechanical Pipettes

Have you ever considered the overall user experience in pipetting? We have.

Sartorius proudly presents Tacta, the new premium mechanical pipette with superb comfort and reliability. Tacta makes pipetting effortless and safe, while producing accurate and reliable results time after time.

At Sartorius, we listened to our customers and combined their views with our extensive experience, and state-of-the-art R&D, to create the exceptionally balanced design of Tacta. Tacta is designed to feel great in your hand, and it is easy and comfortable to use. All materials used for Tacta have been carefully selected, and each component is designed to meet the highest standards.





## The Tacta Experience - Perfectly Balanced

- Comfortable to hold due to the ergonomically designed handle
- Low pipetting and tip ejection forces that reduce the risk of Work Related Upper Limb Disorder (WRULD)
- The unique Sartorius Optilock system provides flexibility for volume adjustment and locking
- Controlled and smooth tip ejection with the new Sartorius Optiject technology
- Calibration adjustment to provide accurate results for various liquid types
- Easy to clean, with only three parts to disassemble
- Reliable results, even over lengthy pipetting periods

# Feel the Comfort

## Superior Ergonomics

Tacta rests lightly in your hand thanks to its ergonomically designed handle and finger hook.

Tacta is effortless to use and exceptionally comfortable to hold. Tacta eases your workload and protects you from strain, even when you pipette for extended periods of your working day. The unique handle and finger hook design lets the pipette rest lightly in your hand, with no need to grip the handle tightly.

## Low Pipetting, Tip Loading and Ejection Forces



Optiload

The forces needed for a pipetting cycle, from tip attachment to pipetting and finally tip ejection, are exceptionally low with Tacta, protecting you from WRULD.

The Optiload feature, with spring-loaded tip cones in both the single and multichannel models, secures tip loading with perfect sealing and minimal

force. Consequently, Optiload also ensures light tip ejection to protect your hand.

The light tip ejection is further enhanced by Sartorius Optiject, the soft, levered tip ejection feature that lets the tip detach in a smooth, controlled way with little force.

## Large, Clear Display



The volume is easy to read even when the pipette is angled.

The large, easy-to-read display helps you to see all four digits of the set volume from a distance without straining your eyes. The volume is easy to read even when the pipette is angled – eliminating the need to turn your head into an uncomfortable position.





## Calibration Adjustment

Tacta is very easy to adjust, for various types of liquids, using a simple adjustment key. The adjustment scale shows the degree of adjustment. By noting this value for a specific liquid, you can return to that setting any time. Easy calibration adjustment ensures accurate results with a variety of liquids.



Calibration Adjustment

## Volume Adjustment

Use the Sartorius Optilock feature to prevent accidental volume changes during pipetting. Either hold down the volume lock button with one hand and adjust the volume with the other, then release to lock, or with just one hand, slide the volume lock up, adjust the volume and slide down to lock again. The choice is in your hands.



Optilock



Volume setting



# Safe from Contamination

## Safe-Cone Filters



Safe-Cone Filters protect the pipette from contamination.


The exchangeable Safe-Cone Filters, used in pipette tip cones, act as barriers to prevent sample aerosols and fluids from contaminating the internal components of the pipette.

The unique Safe-Cone Filter ejector enables the removal of used filters without human contact with contaminated filters.

Safe-Cone Filters are available for all Tacta models greater than 10  $\mu$ l and offer a cost-effective method to reduce contamination. Filters should be changed regularly, and always after over-aspiration.


## Easy Cleaning

Tacta is the quickest and easiest pipette on the market to clean with only three parts to be disassembled for cleaning and no tools are needed for disassembly. Tacta can also be steam-sterilized or autoclaved as-is, without disassembly. It has also high UV and chemical resistance.



The clever design of the Tacta requires only three parts to be disassembled for cleaning, making the Tacta's disassembly time significantly shorter than its competitors'. This makes pipette cleaning an easy and convenient task for anyone and simplifies routine maintenance. Regular cleaning will ensure reliable results and prolong the life of the pipette.


The test was carried out using 1000  $\mu$ l pipettes.



# Achieve Reliable Results

## Proven Precision

The following graph illustrates how the precision of pipetting results is impacted over time when using various pipettes. Tacta, which is designed down to the last detail to create a comfortable pipetting experience, produced consistently precise results during a three-hour test period.



Precision of results tested over a three-hour period using four pipettes from various manufacturers. The test was performed by an independent laboratory in accordance with the ISO 8655 standard using 100 $\mu$ l pipettes with 1000 $\mu$ l tips from the corresponding manufacturers.

# Ordering Information

## Tacta

Order Code	Channels	Volume Range (µl)	Safe-Cone Filters		Optifit Tips* (µl)	SafetySpace Filter Tips* (µl)
			Standard	Plus		
ES81952-103	1	● 0.1-3	–	–	10	10
ES81952-110	1	● 0.5-10	–	–	10	10
ES81952-120	1	● 2-20	721014	–	200	20
ES81952-1100	1	● 10-100	721008	721018	200, 350	120
ES81952-1200	1	● 20-200	721007	721017	200, 350	200
ES81952-11M	1	● 100-1000	721006	721016	1000	1000
ES81952-15M	1	● 500-5000	721005	721015	5000	–
ES81952-110M	1	● 1000-10000	721005	721015	10 000	–
ES81952-810	8	● 0.5-10	–	–	10	10
ES81952-8100	8	● 5-100	721008	721018	200, 350	120
ES81952-8300	8	● 30-300	721007	721017	350	300
ES81952-1210	12	● 0.5-10	–	–	10	10
ES81952-12100	12	● 5-100	721008	721018	200, 350	120
ES81952-12300	12	● 30-300	721007	721017	350	300

\* Low Retention Tips are available in volumes up to 1200µl.



Fully autoclavable

Azer Scientific Inc  
701 Hemlock Road  
Morgantown, PA 19543

Phone: (610) 524-5810 / 1-877-770-2937  
Fax: (610) 901-3046

[www.azersci.com](http://www.azersci.com)