

Rapid-Freeze™

Rapid Freezing Spray



Shown with optional Rapid Zap™ nozzle attachment.

Instantly Cool or Freeze Routine Histology Samples.

Features and Benefits:

Rapid Freeze™ is a non-flammable, moisture-free spray used to instantly cool or freeze samples for routine histology.

Rapid Freeze is able to bring samples to minus 60 degrees Celsius in under three seconds.

Rapid Freeze contains no CFC's or HFC's so it will not deplete stratospheric ozone.

Proprietary *Rapid Zap*™ device is available. It provides a more precise distribution of freeze spray versus standard applicator.

Catalog Number:

RAP-FRZ
RAP-ZAP

Item:

10 oz. (280 grams)
1-Applicator

Packaging:

10 oz. aerosol can
12 per case

For Additional Technical Information call:

EVER Scientific™

Marketing professionals assisting your distribution needs.

337 Gordon Drive • Exton, PA 19341 • P(610)524-5810 • F(610)524-5814 • info@everscientific.com

✓ **Instant Chilling**

✓ **Ozone Safe**

✓ **Precise Distribution**