

● 24 Hour Urine Collection Container

This container is ergonomically designed for refrigeration. It is an HDPE container with a PP cap, which is designed to prevent leakage, further reducing potential contact with urine samples. It also features a precise scale so you may conduct volume measurements quickly and easily. Furthermore, the container is an amber color to protect light-sensitive specimens.

- HDPE Containers with Polypropylene Cap
- 3L Graduated Volume (100mL increments)
- Sized for Easy Refrigeration
- Ergonomic Design



Catalog #

ES3000

Description

24 - Hour Urine Container, Graduated

Packaging

40/case