



AZER SCIENTIFIC

189 Twin County Road Morgantown, PA 19543  
P: 610-524-5810 F: 610-901-3046  
info@azersci.com



## 3 Dimensional Rocker

- *Optimized Speed & Tilt for Gel Blotting, Washing & Staining*
- *Tri-Directional Motion for Thorough, Gentle Mixing*
- *Stacking Platforms Available*
- *Incubator & Cold Room Safe*

### Product Information

Azer Scientific's platform rockers provide the perfect 3-dimensional motion for all blot related applications. Their tilt angle and speed have been optimally set for gentle, but thorough mixing in gel trays, boxes and other flat vessels.

Manufactured in a standard and "mini" version, each rocker comes complete with autoclavable flat mats that are ideal for blot washes and gel staining/destaining. Both models have been designed to conserve bench space. The "mini" rocker is supplied with a 10.5 x 7.5 inch platform while the standard rocker's 12 x 12 inch platform expands the useable work surface by 80%.

### Product Specifications

Max Speed	12rpm
Tilt Angle	+/-8 degrees
Max. Capacity	2.0kg / 4.4 lbs.
Operating Temp	+4°C to +65°C
Weight	2.0kg / 4.4 lbs
Dimension	Mini = 10.5 x 7.5 x 8in    Standard = 12 x 12 x 8.3in
Electrical	115V~, 60 Hz



### Ordering Information

Catalog Number	Item Description
ES3357	Mini 3D Rocker, Fixed Speed, 10" x 7.5" platform w/flat mat
ES3377	Stacking Platform, 10.5" x 7.5" w/flat mat (2.75" separation)
ES3379	Flat Mat, Small, 10.5" x 7.5"
ES3380	Dimpled Mat, Small, 10.5" x 7.5"

Catalog Number	Item Description
ES3358	Standard 3D Rocker, Fixed Speed, 12" x 12" platform w/flat mat
ES3372	Stacking Platform, 12" x 12" w/flat mat (2.75" separation)
ES3374	Flat Mat, Large, 12" x 12"
ES3375	Dimpled Mat, Large, 12" x 12"