



AZER SCIENTIFIC

189 Twin County Road Morgantown, PA 19543
P: 610-524-5810 F: 610-901-3046
info@azersci.com



3 Dimensional Shaker

- *Nutating Motion for Thorough, Uniform Mixing*
- *Dimpled Mat, Perfect for Blood and Centrifuge Tubes*
- *Stacking Platforms Available*
- *Incubator & Cold Room Safe*

Product Information

Azer Scientific's platform shakers provide a 3-dimensional "nutating" motion for moderately aggressive mixing. Two platform sizes and two rubber mat configurations allow for thorough, foam-free mixing in blood tubes and centrifuge tubes, as well as in a variety of flat vessels.

Manufactured in a standard and "mini" version, each shaker comes complete with autoclavable dimpled mats for tube applications. Optional flat mats are available for mixing other vessel types. Both models have been designed to conserve bench space. The "mini" shaker is supplied with a 10.5 x 7.5 inch platform while the standard shaker's 12 x 12 inch platform expands the useable work surface by 80%. Stacking platforms are also available.

Product Specifications

Max Speed	24rpm
Tilt Angle	+/- 20 degrees
Max. Capacity	2.0kg / 4.4 lbs.
Operating Temp	+4°C to +65°C
Weight	2.0kg / 4.4 lbs
Dimension	Mini = 10.5 x 7.5 x 8in Standard = 12 x 12 x 8.3in
Electrical	115V~, 60 Hz



Ordering Information

Catalog Number	Item Description
ES3359	Mini 3D Shaker, Fixed Speed, 10" x 7.5" platform w/flat mat
ES3378	Stacking Platform, 10.5" x 7.5" w/dimpled mat (2.75" separation)
ES3379	Flat Mat, Small, 10.5" x 7.5"
ES3380	Dimpled Mat, Small, 10.5" x 7.5"

Catalog Number	Item Description
ES3361	Standard 3D Shaker, Fixed Speed, 12" x 12" platform w/flat mat
ES3373	Stacking Platform, 12" x 12" w/dimpled mat (2.75" separation)
ES3374	Flat Mat, Large, 12" x 12"
ES3375	Dimpled Mat, Large, 12" x 12"