

# Disposable Petri Dishes

*Superior Clarity, Optimum Flatness, Consistent Quality*

Azer Scientific Petri dishes are manufactured from virgin crystal grade polystyrene for optimum flatness and high optical clarity. Petri dishes are available in a variety of sizes and styles to accommodate both specialized and routine needs. The stacking ring features, on some models, facilitates positive stacking while the slip-ring feature, on other models, facilitates lateral or automated handling. Dishes are sterilized by gamma irradiation and are packaged in linear tear, medical style bags.



**Value**

**Reliability**

**Performance**

**Quality**

**Azer Scientific**

189 Twin County Rd.  
Morgantown, PA 19543

P: 610-524-5810 F: 610-901-3046

[www.azersci.com](http://www.azersci.com)  
[info@azersci.com](mailto:info@azersci.com)

## Standard Petri Dishes

### 60 x 15mm Petri Dish - For Anaerobic Use (Non-Vented)

Cat #	Description	QTY/CS	UM	CS/Pallet	CS Weight
ES3510	60 x 15mm, Sterile, Slipable, No Breather Log	500	25 sleeves, 20 units each	72	11
ES3508	60 x 15mm, Sterile, Semi Stackable, No Breather Log	500	25 sleeves, 20 units each	72	11
ES3506	60 x 15mm, Sterile, Fully Stackable, No Breather Log	500	25 sleeve, 20 units each	72	11

### 60 x 15mm Petri Dish - For Aerobic Use (Vented)

ES3512	60 x 15mm, Sterile, Slipable, w/ Breather Log	500	25 sleeves, 20 units each	72	11
ES3514	60 x 15mm, Sterile, Semi Stackable, w/ Breather Log	500	25 sleeves, 20 units each	72	11
ES3515	60 x 15mm, Sterile, Semi Stackable, w/ Breather Log, Small Lid Sidewall	500	25 sleeves, 20 units each	72	11
ES3505	60 x 15mm, Sterile, Fully Stackable, w/ Breather Log	500	25 sleeves, 20 units each	72	11

### 100 x 15mm, Petri Dish

ES3501	100 x 15mm, Sterile, Mono Plate (Undivided Bottom), Slipable	500	25 sleeves, 20 units each	34	24
ES3561	100 x 15mm, Sterile, Mono Plate (Undivided Bottom), Semi Stackable	500	25 sleeves, 20 units each	34	23
ES3571	100 x 15mm, Sterile, Mono Plate (Undivided Bottom), Fully Stackable	600	25 sleeves, 24 units each	34	27
ES3516	100 x 15mm, Sterile, Mono Plate (Undivided Bottom), Fully Stackable (New Style)	500	25 sleeves, 20 units each	34	24

### 150 x 15mm, Large Petri Dish

ES3525	150 x 15mm, Sterile, Mono Plate (Undivided Bottom)	500	25 sleeves, 20 units each	60	14
--------	--	-----	---------------------------	----	----

## Partitioned Petri Dishes

### 100 x 15mm, Petri Dish - Slipable

ES3502	100 x 15mm, Sterile, I-Plate (2-section), Slipable	500	25 sleeves, 20 units each	34	24
ES3503	100 x 15mm, Sterile, Y-Plate (3-section), Slipable	500	25 sleeves, 20 units each	34	25
ES3504	100 x 15mm, Sterile, X-Plate (4-section), Slipable	500	25 sleeves, 20 units each	34	25

### 100x 15mm, Petri Dish - Semi Stackable

ES3562	100 x 15mm, Sterile, I-Plate (2-section), Semi Stackable	500	25 sleeves, 20 units each	34	23
--------	--	-----	---------------------------	----	----

### 100x 15mm, Petri Dish - Fully Stackable

ES3572	100 x 15mm, Sterile, I-Plate (2-section), Fully Stackable	600	25 sleeves, 24 units each	34	27
ES3573	100 x 15mm, Sterile, Y-Plate (3-section), Fully Stackable	600	25 sleeves, 24 units each	34	28
ES3574	100 x 15mm, Sterile, X-Plate (4-section), Fully Stackable	600	25 sleeves, 24 units each	34	28

### 100x 15mm, Petri Dish - Fully Stackable (New Style)

ES3517	100 x 15mm, Sterile, I-Plate (2-section), Fully Stackable	500	25 sleeves, 20 units each	34	24
ES3518	100 x 15mm, Sterile, Y-Plate (3-section), Fully Stackable	500	25 sleeves, 20 units each	34	25
ES3519	100 x 15mm, Sterile, X-Plate (4-section), Fully Stackable	500	25 sleeves, 20 units each	34	25

## Specialty Petri Dishes

### Contact Plate - 65 x 15mm

Cat #	Description	QTY/CS	UM	CS/Pallet	CS Weight
ES3550	65 x 15mm, Contact Plate, Sterile, Flat Bottom, Classic Style	500	25 sleeves, 20 units each	64	13
ES3556	65 x 15mm, Contact Plate, Sterile, Convex Bottom, Outside Grid	500	25 sleeves, 20 units each	64	13
ES3586	65 x 15mm, Contact Plate, Sterile, Convex Bottom, Center Gated, Inside Grid	500	25 sleeves, 20 units each	64	13
ES3577	65 x 15mm, Contact Plate, Sterile, Convex Bottom, Center Grated, Outside Grid	700	25 sleeves, 28 units each	64	17

### Locking Contact Plate - 65 x 15mm

ES3585	65 x 15mm, Contact Plate, Sterile, Classic Style, Convex Bottom, Friction Lock Design	700	25 sleeves, 28 units each	64	17
ES3598	65 x 15mm, Contact Plate, Sterile, Classic Style, Convex Bottom, Rotation Lock Design	700	25 sleeves, 28 units each	64	17

### Locking Contact Plate - 100 x 10mm

ES3511	100 x 10mm, Sterile, Mono Plate (Undivided Bottom)	750	25 sleeves, 30 units each	34	29
--------	--	-----	---------------------------	----	----

### Petri Dish for Anaerobic Use with CO<sub>2</sub> Tablet

ES3596	75 x 50mm, Sterile, Rectangular, With Pill Pocket for CO <sub>2</sub> Tablet	500	25 sleeves, 20 units each	64	14
ES3509	100 x 15mm, Sterile, Mono Plate (Undivided Bottom), Slipable, With Pill Pocket for CO <sub>2</sub> Tablet	500	25 sleeves, 20 units each	34	24
ES3552	100 x 15mm, Sterile, I-Plate (2-section), Slipable, With Pill Pocket for CO <sub>2</sub> Tablet	500	25 sleeves, 20 units each	34	24

## Newly Added Products

These dishes are ideal for phage typing, susceptibility testing, plate counts, and probe assays. Each 13mm grid is marked numerically in one direction and alphabetically in the other. This dish has four venting ribs in the underside of the lid to prevent condensation build-up.



<u>Cat #</u>	<u>Description</u>	<u>UM</u>
ES3308	Petri Dish, Square w/ grid, Sterile, 100 x 100 x 15, PS, Vented	10/ sleeve, 500/cs

These dishes accommodate deeper fills for longer culture periods. Used for fungal cultures and plant propagation. This dish comes vented and non-vented.



<u>Cat #</u>	<u>Description</u>	<u>UM</u>
ES3306	Petri Dish, Sterile, 125ml, 100 x 25, PS, Vented	20 per Sleeve, 500/Case
ES3307	Petri Dish, Sterile, 125ml, 100 x 25, PS, Non-Vented, Stackable	20 per Sleeve, 500/Case

These 50 x 9mm sterile Petri dishes with absorbent pads are ideal for culturing micro-organisms on either agar or broth based media. Designed to accommodate 47mm diameter membrane filters. These dishes are stackable and a frosted area on top of the lid permits labeling.



<u>Cat #</u>	<u>Description</u>	<u>UM</u>
ES3309	Petri Dish w/ Absorbent Pad, Sterile, 50 x 9, PS, Stackable	20 per Sleeve, 500/Case