



AZER SCIENTIFIC

189 Twin County Road—Morgantown, PA 19543
P: (610)-524-5810 F: (610)-901-3046
www.azersci.com

● Tissue and Biopsy Processing Cassettes



Tissue / Biopsy Two Piece Cassettes

Suitable for most histological applications, these two piece cassettes are molded from a special high density polymer that is totally resistant to the chemical action of solvents. Lids can be snapped on and off as often as necessary and they always relock securely without danger of specimen loss. Large labeling areas are located on three sides of the cassettes.

Lids and bases packaged separately in boxes of 500.



Cat #	Description	Qty/Cs
ES3013-B	Embedding Cassette, Blue	500
ES3013-GY	Embedding Cassette, Gray	500
ES3013-GN	Embedding Cassette, Green	500
ES3013-L	Embedding Cassette, Lilac	500
ES3013-O	Embedding Cassette, Orange	500
ES3013-P	Embedding Cassette, Pink	500
ES3013-T	Embedding Cassette, Tan	500
ES3013-W	Embedding Cassette, White	500
ES3013-Y	Embedding Cassette, Yellow	500
ES3014-B	Biopsy Cassette, Blue	500
ES3014-G	Biopsy Cassette, Green	500
ES3014-W	Biopsy Cassette, White	500
ES3014-Y	Biopsy Cassette, Yellow	500

Tissue / Biopsy Cassettes in FastLoad Sleeves

FastLoad Tissue Cassette

Our transparent sleeve allows viewing of cassettes in order to confirm there is no jam in the sleeve during the printing process. A smooth inside finish ensures that cassettes will never jam in the sleeve while the cassette labeler is in use.

These embedding cassettes are designed with efficient flow-through slots to maximize fluid exchange and ensure proper drainage. Lids can be opened and closed as often as necessary and they always relock securely without danger of specimen loss. When ready for the embedding process, covers can conveniently be snapped off and discarded.

Each case contains 10 sleeves of 75 cassettes.



Cat #	Description	Qty/Cs
ES3015-W	Tissue Cassette, FastLoad, White	750
ES3015-P	Tissue Cassette, FastLoad, Pink	750
ES3015-G	Tissue Cassette, FastLoad, Green	750
ES3015-Y	Tissue Cassette, FastLoad, Yellow	750
ES3015-B	Tissue Cassette, FastLoad, Blue	750
ES3015-PH	Tissue Cassette, FastLoad, Peach	750
ES3015-T	Tissue Cassette, FastLoad, Tan	750


Cat #	Description	Qty/Cs
ES3015-GY	Tissue Cassette, FastLoad, Gray	750
ES3015-L	Tissue Cassette, FastLoad, Lilac	750
ES3015-O	Tissue Cassette, FastLoad, Orange	750
ES3015-A	Tissue Cassette, FastLoad, Aqua	750
ES3015-FP	Tissue Cassette, FastLoad, Fluo Pink	750
ES3015-FG	Tissue Cassette, FastLoad, Fluo Green	750
ES3015-FY	Tissue Cassette, FastLoad, Flou Yellow	750

FastLoad Biopsy Cassette

Suitable for hoppers accepting plastic sleeves, these cassettes will load in cassette labeling instruments in one simple operation. Just load the 75-cassette FastLoad sleeve in the hopper and you are ready for printing.

These biopsy cassettes are molded from a special high density polymer that is totally resistant to the chemical action of histological solvents. Lids can be opened and closed as often as necessary and they always relock securely without danger of specimen loss. Large labeling areas are located on three sides of the cassettes.

Each case contains 10 sleeves of 75 cassettes.



Cat #	Description	Qty/Cs
ES3016-W	Biopsy Cassette, FastLoad, White	750
ES3016-P	Biopsy Cassette, FastLoad, Pink	750
ES3016-G	Biopsy Cassette, FastLoad, Green	750
ES3016-Y	Biopsy Cassette, FastLoad, Yellow	750
ES3016-B	Biopsy Cassette, FastLoad, Blue	750
ES3016-PH	Biopsy Cassette, FastLoad, Peach	750
ES3016-T	Biopsy Cassette, FastLoad, Tan	750

Cat #	Description	Qty/Cs
ES3016-GY	Biopsy Cassette, FastLoad, Gray	750
ES3016-L	Biopsy Cassette, FastLoad, Lilac	750
ES3016-O	Biopsy Cassette, FastLoad, Orange	750
ES3016-A	Biopsy Cassette, FastLoad, Aqua	750
ES3016-FP	Biopsy Cassette, FastLoad, Fluo Pink	750
ES3016-FG	Biopsy Cassette, FastLoad, Fluo Green	750
ES3016-FY	Biopsy Cassette, FastLoad, Fluo Yellow	750

FastLoad Tissue Cassette

Suitable for hoppers accepting plastic sleeves, these cassettes will load in cassette labeling instruments in one simple operation. Just load the 75-cassette FastLoad sleeve in the hopper and you are ready for printing.

These mesh cassettes offer 0.38mm square openings. To allow for increased fluid exchange and faster sinking in liquids, these cassettes have large anterior slots. No biopsy pads are necessary. Each cassette also features a recessed cover and a large square 27mm compartment, perfect even for needle biopsies.

Each case contains 10 sleeves of 75 cassettes.



Cat #	Description	Qty/Cs
ES3017-L	Mesh Cassette, FastLoad, Lilac	750
ES3017-O	Mesh Cassette, FastLoad, Orange	750
ES3017-P	Mesh Cassette, FastLoad, Pink	750
ES3017-PH	Mesh Cassette, FastLoad, Peach	750
ES3017-T	Mesh Cassette, FastLoad, Tan	750
ES3017-W	Mesh Cassette, FastLoad, White	750
ES3017-Y	Mesh Cassette, FastLoad, Yellow	750

Cat #	Description	Qty/Cs
ES3017-A	Mesh Cassette, FastLoad, Aqua	750
ES3017-B	Mesh Cassette, FastLoad, Blue	750
ES3017-FG	Mesh Cassette, FastLoad, Fluo Green	750
ES3017-FP	Mesh Cassette, FastLoad, Fluo Pink	750
ES3017-FY	Mesh Cassette, FastLoad, Flou Yellow	750
ES3017-G	Mesh Cassette, FastLoad, Green	750
ES3017-GY	Mesh Cassette, FastLoad, Gray	750

Tissue / Biopsy Hinged Cassettes

Hinged Tissue Cassette

These cassettes are more compact, easier to use and more efficient than ever. Each cassette incorporates a unique recessed cover, a great space saving feature allowing more cassettes to be stacked in automatic labeling machines and in storage cabinets.

Molded from a special high density polymer, these cassettes keep specimens safely submerged in liquid and are totally resistant to the chemical action of histological solvents. These cassettes ensure efficient fluid exchange and drainage thanks to 114 openings each measuring 1 x 5 mm. The one-piece integral lid eliminated the need for separate steel lids. They can be opened and closed as often as necessary and they always relock securely without danger of specimen loss. Large labeling areas for easy identification.

Each case contains 3 dispenser boxes of 500 cassettes.



Cat #	Description	Qty/Cs
ES3336-L	Tissue Cassette, Assembled, Lilac	1500
ES3336-O	Tissue Cassette, Assembled, Orange	1500
ES3336-P	Tissue Cassette, Assembled, Pink,	1500
ES3336-PH	Tissue Cassette, Assembled, Peach	1500
ES3336-T	Tissue Cassette, Assembled, Tan	1500
ES3336-W	Tissue Cassette, Assembled, White	1500
ES3336-Y	Tissue Cassette, Assembled, Yellow	1500

Cat #	Description	Qty/Cs
ES3336-A	Tissue Cassette, Assembled, Aqua	1500
ES3336-B	Tissue Cassette, Assembled, Blue	1500
ES3336-FG	Tissue Cassette, Assembled, Fluo Green	1500
ES3336-FP	Tissue Cassette, Assembled, Fluo Pink	1500
ES3336-FY	Tissue Cassette, Assembled, Fluo Yellow	1500
ES3336-G	Tissue Cassette, Assembled, Green	1500
ES3336-GY	Tissue Cassette, Assembled, Gray	1500

Hinged Biopsy Cassette

Similar to the ES3017 cassette series, but designed specifically to hold biopsy specimens during the embedding process, as well as in a storage cabinet. Anterior writing area is at a 45° angle to make the cassette more suitable to be used with certain types of cassette labeling instruments.

Each case contains 3 dispenser boxes of 500 cassettes.



Cat #	Description	Qty/Cs	Cat #	Description	Qty/Cs
ES3406-L	Biopsy Cassette, Assembled, Lilac	1500	ES3406-A	Biopsy Cassette, Assembled, Aqua	1500
ES3406-O	Biopsy Cassette, Assembled, Orange	1500	ES3406-B	Biopsy Cassette, Assembled, Blue	1500
ES3406-P	Biopsy Cassette, Assembled, Pink	1500	ES3406-FG	Biopsy Cassette, Assembled, Fluo Green	1500
ES3406-PH	Biopsy Cassette, Assembled, Peach	1500	ES3406-FP	Biopsy Cassette, Assembled, Fluo Pink	1500
ES3406-T	Biopsy Cassette, Assembled, Tan	1500	ES3406-FY	Biopsy Cassette, Assembled, Fluo Yellow	1500
ES3406-W	Biopsy Cassette, Assembled, White	1500	ES3406-G	Biopsy Cassette, Assembled, Green	1500
ES3406-Y	Biopsy Cassette, Assembled, Yellow	1500	ES3406-GY	Biopsy Cassette, Assembled, Gray	1500

Hinged Mesh Cassette

These mesh cassettes offer 1676 square openings (0.38 mm) allowing for a greatly improved fluid exchange without having to use biopsy pads. Large anterior and posterior slots in both the cassette and cover ensure that the cassette will sink rapidly. A large square compartment measuring 27mm is perfect even for needle biopsies.

Molded from a special high density polymer, these cassettes keep specimens safely submerged in liquid and are totally resistant to the chemical action of histological solvents. These cassettes ensure efficient fluid exchange and drainage. Anterior writing area is at a 45° angle to make the cassette more suitable to be used with certain types of cassette labeling instruments.

Each case contains 4 dispenser boxes of 250 cassettes.



Cat #	Description	Qty/Cs
ES3335-L	Mesh Cassette, Hinged, Lilac	1000
ES3335-O	Mesh Cassette, Hinged, Orange	1000
ES3335-P	Mesh Cassette, Hinged, Pink	1000
ES3335-PH	Mesh Cassette, Hinged, Peach	1000
ES3335-T	Mesh Cassette, Hinged, Tan	1000
ES3335-W	Mesh Cassette, Hinged, White	1000
ES3335-Y	Mesh Cassette, Hinged, Yellow	1000

Cat #	Description	Qty/Cs
ES3335-A	Mesh Cassette, Hinged, Aqua	1000
ES3335-B	Mesh Cassette, Hinged, Blue	1000
ES3335-FG	Mesh Cassette, Hinged, Fluo Green	1000
ES3335-FP	Mesh Cassette, Hinged, Fluo Pink	1000
ES3335-FY	Mesh Cassette, Hinged, Flou Yellow	1000
ES3335-G	Mesh Cassette, Hinged, Green	1000
ES3335-GY	Mesh Cassette, Hinged, Gray	1000

Hinged Tissue & Biopsy Cassette

These disposable plastic cassettes hold tissue specimens during processing and embedding, as well as in storage. Molded from a special high-density polymer, they keep specimens safely submerged in liquid and are totally resistant to the chemical action of histological solvents. The efficient flow-through slots maximize fluid exchange and ensure proper drainage.



This new model differs by the special hinge holding the base and cover together. This hinge allows cassettes to be opened and closed as often as necessary. The cover can be removed and re-inserted easily without danger of specimen loss. Three frosted writing surfaces are available, including the anterior one at a 45° angle.



Each case contains 3 dispenser boxes of 500 cassettes.

Hinged Tissue Cassettes

Hinged Biopsy Cassettes

Cat #	Description	Qty/Cs	Cat #	Description	Qty/Cs
ES3018-A	Tissue Cassette, Hinged, Aqua	1500	ES3407-A	Biopsy Cassette, Hinged, Aqua	1500
ES3018-B	Tissue Cassette, Hinged, Blue	1500	ES3407-B	Biopsy Cassette, Hinged, Blue	1500
ES3018-FG	Tissue Cassette, Hinged, Fluo Green	1500	ES3407-FG	Biopsy Cassette, Hinged, Fluo Green	1500
ES3018-FP	Tissue Cassette, Hinged, Fluo Pink	1500	ES3407-FP	Biopsy Cassette, Hinged, Fluo Pink	1500
ES3018-FY	Tissue Cassette, Hinged, Flou Yellow	1500	ES3407-FY	Biopsy Cassette, Hinged, Flou Yellow	1500
ES3018-G	Tissue Cassette, Hinged, Green	1500	ES3407-G	Biopsy Cassette, Hinged, Green	1500
ES3018-GY	Tissue Cassette, Hinged, Gray	1500	ES3407-GY	Biopsy Cassette, Hinged, Gray	1500
ES3018-L	Tissue Cassette, Hinged, Lilac	1500	ES3407-L	Biopsy Cassette, Hinged, Lilac	1500
ES3018-O	Tissue Cassette, Hinged, Orange	1500	ES3407-O	Biopsy Cassette, Hinged, Orange	1500
ES3018-P	Tissue Cassette, Hinged, Pink	1500	ES3407-P	Biopsy Cassette, Hinged, Pink	1500
ES3018-PH	Tissue Cassette, Hinged, Peach	1500	ES3407-PH	Biopsy Cassette, Hinged, Peach	1500
ES3018-T	Tissue Cassette, Hinged, Tan	1500	ES3407-T	Biopsy Cassette, Hinged, Tan	1500
ES3018-W	Tissue Cassette, Hinged, White	1500	ES3407-W	Biopsy Cassette, Hinged, White	1500
ES3018-Y	Tissue Cassette, Hinged, Yellow	1500	ES3407-Y	Biopsy Cassette, Hinged, Yellow	1500

Great Alternative to the Richard-Allan and Surgipath Cassettes!

