

Azer Scientific's 15 mL and 50 mL sterile centrifuge tubes are manufactured from high quality clarified polypropylene, non-pyrogenic, and endotoxin free. Each size is available in bulk or racked configuration, feature black graduations and a white labeling area, and can withstand temperatures from -45°C to 121°C (autoclavable).

Features and Benefits

- **Manufactured from high quality polypropylene**
- **Manufactured in an ISO 13485 facility**
- **Black printed graduations and white labeling area**
- **Non-pyrogenic, endotoxin free, ensuring the integrity of your samples**



Azer Scientific carries a complete line of lab products, including:

- Small equipment
- Cryogenic vials
- Liquid handling products
- Microscope slides & cover glass

Visit
www.AzerSci.com
for more!

Specifications:

15 mL

Graduation Range: 1 - 14.5 mL
Temp Range: -45°C - 121°C
Max RCF: 10,000xg

50 mL

Graduation Range: 5 - 47.5 mL
Temp Range: -45°C - 121°C
Max RCF: 9,500xg

| Item | Description | Quantity |
|----------|--|----------|
| ES-0151N | 15 mL centrifuge tube, plug cap, sterile, racked | 500/case |
| ES-0152N | 15 mL centrifuge tube, plug cap, sterile, bulk | 500/case |
| ES-9151N | 15 mL centrifuge tube, flat cap, sterile, racked | 500/case |
| ES-9152N | 15 mL centrifuge tube, flat cap, sterile, bulk | 500/case |
| ES-0501N | 50 mL centrifuge tube, plug cap, sterile, racked | 500/case |
| ES-0502N | 50 mL centrifuge tube, plug cap, sterile, bulk | 500/case |
| ES-9501N | 50 mL centrifuge tube, flat cap, sterile, racked | 500/case |
| ES-9502N | 50 mL centrifuge tube, flat cap, sterile, bulk | 500/case |