

Item #	Description
ESCS-102	Control Slide, Acid Fast Bacteria
ESCS-103	Control Slide, Pnuemocystis
ESCS-104	Control Slide, Alcian Blue
ESCS-105	Control Slide, Amoeba
ESCS-106	Control Slide, Amyloid
ESCS-107	Control Slide, Argentaffin
ESCS-108	Control Slide, Austalian Antigen
ESCS-109	Control Slide, Bile
ESCS-110	Control Slide, Calcium
ESCS-111	Control Slide, Elastic Tissue
ESCS-112	Control Slide, Fat
ESCS-113	Control Slide, Fungi-Pas
ESCS-114	Control Slide, Fungi-Gms Silver Strain
ESCS-115	Control Slide, Giemsa
ESCS-116	Control Slide, Glycogen
ESCS-117	Control Slide, Gram Organism
ESCS-118	Control Slide, Iron
ESCS-119	Control Slide, Leprosy
ESCS-120	Control Slide, Luxol Fast Blue
ESCS-121	Control Slide, Massons
ESCS-122	Control Slide, Melanin
ESCS-123	Control Slide, Methyl Green Pyronin
ESCS-124	Control Slide, Muci Carmine
ESCS-125	Control Slide, Periodic Acid
ESCS-126	Control Slide, Ptah
ESCS-127	Control Slide, Reticulum
ESCS-128	Control Slide, Spirochetes
ESCS-129	Control Slide, Trichrome
ESCS-130	Control Slide, Helicobacter

Introducing our Full Line of Histology Control Slides

Azer Scientific's control slides are carefully prepared with paraffin embedded human tissue sections containing specific characteristics that are used for monitoring staining procedures in pathology laboratories. Each package contains 10 EverMark microscope slides. One slide is stained with the proper procedure to demonstrate and validate the expected result.




Photo of Muci Carmine (ESCS124)

AZER SCIENTIFIC
 189 Twin County Road
 Morgantown, PA 19543
 P: 610-524-5810 F: 610-901-3046
info@AzerSci.com