




AZER SCIENTIFIC

● CYTOLOGY SPIN SLIDES

- Compatible with the Shandon™ Centrifuge Systems.
- Designed to make locating the centrifuged smear quicker and easier.
- Each well is placed in the exact position where the centrifuge specimen is deposited.
- Available with single or double rings, with or without adhesive coating.



ESCS-096	Cytology Slide, Shandon™ centrifuge systems. Single Well, Uncoated, per Gross
ESCS-097	Cytology Slide, Shandon™ centrifuge systems. Double Well, Uncoated, per Gross
ESCS-098	Cytology Slide, Shandon™ centrifuge systems. Single Well Positively Charged, per Gross
ESCS-099	Cytology Slide, Shandon™ centrifuge systems. Double Well Positively Charged, per Gross



Azer Scientific

189 Twin County Road
Morgantown, PA 19543

P: 610-524-5810

F: 610-901-3046

info@azersci.com