



AZER SCIENTIFIC

189 Twin County Rd. - Morgantown, PA 19543

P: 610-524-5810 F: 610-901-3046

www.azersci.com



# Drosophila Products

*Azer Scientific manufactures Drosophila Products that offer flexibility and product excellence at the best price. These products are designed specifically for Drosophilae Research.*



Our bottles are reusable yet disposably priced, manufactured with thin walls for flexibility. They are available in transparent autoclavable polypropylene or less transparent non-autoclavable polyethylene.



Bottle openings are designed to accept standard size paper lids, foam plugs, cotton balls, and Buzz Plugs™.



Vials are available in either standard cost saving polystyrene or flexible virtually unbreakable K-resin. Our vials have a flat bottom straight wall design for convenient racking, stacking and filling.



Vials are available in tray or bulk pack. The tray pack 10 x 10 row configuration is great for quick, effortless filling. Bulk pack vials are layered with open-end orientation for convenient handling and smaller case size.

## Bottles

Cat #	Description	Dimensions	CS/Weight	QTY/CS
ES-38115	Round bottom, Polypropylene, 6oz with indent for Buzz Plug	64 diam. X 103 high	21 lbs	500
ES-38354	Square bottom, Polypropylene, 6oz with wide ledge for paper lid	57 x 57 x 103 high	21 lbs	500
ES-38355	Square bottom, Polypropylene, 6oz with indent for Buzz Plug	57 x 57 x 103 high	18 lbs	500
ES-38359	Square bottom, Polyethylene, 6oz, with indent for Buzz Plug	57 x 57 x 103 high	21 lbs	500
ES-38362	Square bottom, Polyethylene, 6oz, with wide ledge for paper lid	57 x 57 x 103 high	18 lbs	500



## Vials

Cat #	Material	Style	Pack	Size (mm)	Weight/CS	QTY/CS
ES-3825	Polystyrene	Shorty	Bulk	24 x 75	9 lbs	500
ES-3826	Polystyrene	Shorty	Tray	24 x 75	11 lbs	500
ES-3815	Polystyrene	Narrow	Bulk	25 x 95	10 lbs	500
ES-3816	Polystyrene	Narrow	Tray	25 x 95	11 lbs	500
ES-3819	Polystyrene	Wide	Bulk	28.5 x 95	12 lbs	500
ES-3820	Polystyrene	Wide	Tray	28.5 x 95	13 lbs	500
ES-3809	K-Resin	Narrow	Bulk	25 x 95	9 lbs	500
ES-3808	K-Resin	Narrow	Tray	25 x 95	10 lbs	500
ES-3811	K-Resin	Wide	Bulk	28.5 x 95	11 lbs	500
ES-3810	K-Resin	Wide	Tray	28.5 x 95	12 lbs	500
ES-3807	Polypropylene	Narrow	Bulk	25 x 95	10 lbs	500
ES-3814	Polypropylene	Narrow	Tray	25 x 95	13 lbs	500
ES-3813	Polypropylene	Wide	Bulk	28.5 x 95	11 lbs	500
ES-3812	Polypropylene	Wide	Tray	28.5 x 95	13 lbs	500



## Petri Dishes

Cat #	Size (mm)	Description	Weight/CS	QTY/CS
ES3508	60 x 15	Non-vented with stacking rim	10 lbs	20/sleeve, 500/cs
ES3506	60 x 15	Non-vented with rounded edge & stacking rim	10 lbs	20/sleeve, 500/cs
ES3510	60 x 15	Non-vented with rounded edge & no stacking rim	10 lbs	20/sleeve, 500/cs
ES3505	60 x 15	Vented with stacking rim	10 lbs	20/sleeve, 500/cs