

## Pre-Fill Formalin Containers

*This is a traditional fixative that is buffered to a neutral pH. It is used in histological applications for fixing and storing tissue specimens. This line of all plastic containers was designed specifically for the collection and transport of various sizes of tissue samples. Each container is half filled with 10% neutral buffered formalin. The most commonly used sizes: 20, 40, 60, 90, and 120ml varieties are packaged with our proprietary Sure-Lok system. These containers are 95 kPa compliant.*

95 kPa Compliant Containers

Sure-Lok Sealing System



We carry 8 place, 12 place and 16 place Prostate Biopsy Kits. For more information visit us on the web, [www.azersci.com](http://www.azersci.com).

Catalog #	Description	UOM
PFNBF-20	20ml Container with 10ml 10% Neutral Buffered Formalin	128 per case
B-PFNBF-20	20ml Container with 10ml 10% Neutral Buffered Formalin - Bulk	784 per case
PFNBF-40	40ml Container with 20ml 10% Neutral Buffered Formalin	50 per case
B-PFNBF-40	40ml Container with 20ml 10% Neutral Buffered Formalin - Bulk	420 per case
PFNBF-60	60ml Container with 30ml 10% Neutral Buffered Formalin	50 per case
B-PFNBF-60	60ml Container with 30ml 10% Neutral Buffered Formalin - Bulk	300 per case
PFNBF-90	90ml Container with 45ml 10% Neutral Buffered Formalin	50 per case
B-PFNBF-90	90ml Container with 45ml 10% Neutral Buffered Formalin - Bulk	240 per case
PFNBF-120	120ml Container with 60ml 10% Neutral Buffered Formalin	50 per case
B-PFNBF-120	120ml Container with 60ml 10% Neutral Buffered Formalin - Bulk	200 per case
PFNBF-180	180ml Container with 90ml 10% Neutral Buffered Formalin	50 per case
PFNBF-240	240ml Container with 120ml 10% Neutral Buffered Formalin	24 per case
PFNBF-360	360ml Container with 180ml 10% Neutral Buffered Formalin	24 per case
PFNBF-480	480ml Container with 240ml 10% Neutral Buffered Formalin	24 per case
PFNBF-1000	1000ml Container with 500ml 10% Neutral Buffered Formalin	12 per case
PFNBF-3840	3840ml Container with 1920ml 10% Neutral Buffered Formalin	4 per case
NBF-4-G	Buffered 10% Formalin Fixative, 1 gallon	4 per case
CUNBF-2.5G	Buffered 10% Formalin Fixative, 2.5 gallon	1 cube
CUNBF-5-G	Buffered 10% Formalin Fixative, 5 gallon	1 cube
NBF-30G	Buffered 10% Formalin Fixative, 30 gallon	1 drum
NBF-55G	Buffered 10% Formalin Fixative, 55 gallon	1 drum



AZER SCIENTIFIC

189 Twin County Road - Morgantown, PA 19543 - P: 610-524-5810 - F: 610-901-3046 - info@azersci.com