



Large Pathology Containers

for secure storage and transportation of large specimens

Our large pathology containers are sturdy, durable and shatter resistant. A screw-top lid creates a secure seal that helps prevent leakage. Manufactured from high-density polyethylene, these containers are a great storage option for your lab. These containers feature a latching closure that keeps lids secure and a carry handle for easy mobility.



Cat. #	Description	Quantity
MM100	Large Pathology Container with screw-top lid, 0.6 gallon	12/case
MM200	Large Pathology Container with screw-top lid, 1.25 gallon	12/case
MM350	Large Pathology Container with screw-top lid, 3.5 gallon	4/case
MM500	Large Pathology Container with screw-top lid , 5 gallon	4/case

