




AZER SCIENTIFIC

189 Twin County Road  
Morgantown, PA 19543  
P: 610-524-5810 F: 610-901-3046  
info@azersci.com  
www.azersci.com

# Pathology Containers

*for Storage of Lab Specimens*

Our Pathology Containers are sturdy, light-weight and shatter resistant. A snap-on lid creates a secure seal that helps prevent leakage. Manufactured from translucent high-density polyethylene, these containers have excellent resistance to acids, bases and alcohols. Temperature range: freezing - boiling.



CAT #	DESCRIPTION	"A" [in]	"B" [in]
ES3700	Pathology Container (non-sterile) with lid- 4 oz, 250/case	3	2-1/4
ES3701	Pathology Container (non-sterile) with lid- 8 oz, 250/case	3-1/4	2-5/8
ES3701-S	Pathology Container (non-sterile) with lid- 8 oz, 250/case	3-7/8	2-1/4
ES3702	Pathology Container (non-sterile) with lid- 16 oz, 100/case	3-3/4	3-3/4
ES3703	Pathology Container (non-sterile) with lid- 32 oz, 100/case	4-3/4	4-3/4
ES3704	Pathology Container (non-sterile) with lid- 64 oz, 50/case	6-1/2	4-3/4
ES3705	Pathology Container (non-sterile) with lid- 84 oz, 25/case	6-1/2	6-1/2
ES3706	Pathology Container (non-sterile) with lid- 165 oz, 25/case	8-1/4	8