



Rapid Freeze™

INSTANTLY COOL OR FREEZE HISTOLOGY SAMPLES

Features & Benefits

- *Rapid Freeze™* is a non-flammable, moisture free spray used to instantly cool or freeze samples for routine histology.
- *Rapid Freeze™* is able to bring sample to minus 60 degrees Celsius in under three seconds.
- *Rapid Freeze™* contains no CFC's or HFC's so it will not deplete stratospheric ozone.
- Proprietary *Rapid Zap™* provides a more precise distribution of freeze spray versus the standard applicator.



Catalog

RAP-FRZ

RAP-ZAP

Description

Rapid Freeze, 10oz aerosol can

Applicator nozzle for Rapid Freeze can

UM

12/case

1 each

Instant Chilling

Ozone Safe

Precise Distribution