

Slide and Tissue Block Storage File



- Economical
- Stackable, Interlocking Construction
- Easy to Assemble
- Box and Drawers Allow Indexing
- Systems have a Polymer-Coated Surface that Helps to Resist Moisture
- Dimensions = 12" wide x 11 1/2" long x 4" tall

<u>Catalog #</u>	<u>Description</u>	<u>UOM</u>
ES040	4 - Drawer Slide File, accommodates 2000 slides	10 per case
ES080	8 - Drawer Tissue Block File, Accommodates 500 blocks	10 per case

