



Vortex Mixer

- *Powerful Motor for INSTANT Vortexing*
- *Q-Drive™ Dynamic Balance System*
- *Touch or Continuous Operation*
- *Variable Speed Control from 200 to 3200rpm*

Product Information

Azer Scientific's vortex mixer sets the new standard in vortex mixing. With advanced Q-Drive™ technology, our vortex mixer provides smooth instant vortexing of even the largest sample sizes, including full 50ml tubes. The unique counter balance system creates maximum vortexing action, while minimizing noise and excessive vibration. This results in quieter, more efficient and longer lasting operation.

Optional combination head available. Holds 66 tubes (38 x 1.5ml and 28 x 0.5ml) or one microplate.

Product Specifications

Max Speed	200 to 3200rpm
Operation Mode	Touch or continuous
Operating Temp	+4°C to +65°C
Weight	3.8kg / 8.4 lbs
Dimension	5 x 6.3 x 6.7" (13 x 16 x 17cm)
Electrical	115V~, 60 Hz

“No-Walk” Suction Feet

Cold Room/Incubator Compatibility

Spill-Proof Electronic Controls

Ordering Information

Catalog Number	Item Description
ES3362	Vortex Mixer
ES3370	Optional 3" flat head
ES3371	Optional combination head



189 Twin County Road
Morgantown, PA 19543
P: 610-524-5810 F: 610-901-3046
info@azersci.com