

Newly Added Products

These dishes are ideal for phage typing, susceptibility testing, plate counts, and probe assays. Each 13mm grid is marked numerically in one direction and alphabetically in the other. This dish has four venting ribs into the underside of the lid to prevent condensation build-up.

Cat #

Description

UM

ES3308 Petri Dish, Square w/ grid, Sterile, 100 x 15, PS, Vented

10/ sleeve, 500/cs



These dishes accommodate deeper fills for longer culture periods. Used for fungal cultures and plant propagation. These dishes come vented and non-vented.

Cat #

Description

UM

ES3306 Petri Dish, Sterile, 125ml, 100 x 25, PS, Vented

20 per Sleeve, 500/Case

ES3307 Petri Dish, Sterile, 125ml, 100 x 25, PS, Non-Vented, Stackable

20 per Sleeve, 500/Case



These 50 x 9mm sterile Petri dishes with absorbent pads are ideal for culturing micro-organisms on either agar or broth based media. Designed to accommodate 47mm diameter membrane filters. These dishes are stackable and a frosted area on top of the lid permits labeling.

Cat #

Description

UM

ES3309 Petri Dish w/ Absorbent Pad, Sterile, 50 x 9, PS, Stackable

20 per Sleeve, 500/Case

