

Certificate of Analysis

Product Name: Pre-Filled Neutral Buffered Formalin, 240ml

Description: 10% Neutral Buffered Formalin, Prefilled Container, 240ml, Half filled, 24/cs

Product No.: PFNBF-240

Lot No.: 26466

Date of Manufacture: 10/27/2021

Expiration Date: 2023/10/01

Materials:

	<u>Part No.</u>	<u>Lot No.</u>	<u>Exp. Date</u>
Chemical:	NBF-275F	1J010	2023/10/01

QC Testing:

LOT NUMBER: 1J010			
PROPERTY	SPECIFICATION	RESULTS	PASS / FAIL
pH	6.90-7.10 @ 25°C	7.03 @ 25°C	PASS
Formaldehyde Concentration	9.5-10.5% (v/v)	9.70 % (v/v)	PASS
Appearance	Clear with No Particulate	Clear	PASS
Labeling	Correct Lot Number and Exp. Date	Correct	PASS

Country of Origin: United States of America

Latex Free: This product contains no latex

Created By:

Fernando Fana (QAA)