

October 26, 2021

CERTIFICATE OF ANALYSIS

Product Name: 10% Neutral Buffered Formalin, Gallon
Description: 10% Neutral Buffered Formalin, Gallon, 4 x 1 gallon/case
Product No.: NBF-4-G
Lot No.: 26698
Date of Manufacture: 2021-10-21
Expiration Date: 2023-10-01

Materials:

	<u>Part No.</u>	<u>Lot No.</u>	<u>Exp. Date</u>
Chemical:	NBF-275G	1J005	2023-10-01
	NBF-275G	1J006	2023-10-01

QC Testing:

LOT NUMBER: 1J005			
PROPERTY	SPECIFICATION	RESULTS	PASS / FAIL
pH	6.90-7.10 @ 25°C	7.02 @ 25°C	PASS
Formaldehyde Concentration	9.5-10.5% (v/v)	10.31% (v/v)	PASS
Appearance	Clear with No Particulate	Clear	PASS
Labeling	Correct Lot Number and Exp. Date	Correct	PASS
LOT NUMBER: 1J006			
PROPERTY	SPECIFICATION	RESULTS	PASS / FAIL
pH	6.90-7.10 @ 25°C	7.06 @ 25°C	PASS
Formaldehyde Concentration	9.5-10.5% (v/v)	10.31% (v/v)	PASS
Appearance	Clear with No Particulate	Clear	PASS
Labeling	Correct Lot Number and Exp. Date	Correct	PASS

Country of Origin: United States of America

Latex Free: This product contains no latex

 Created By:
 James Osborn
 Quality Associate