

October 27, 2021

### CERTIFICATE OF ANALYSIS

**Product Name:** 10% Neutral Buffered Formalin, 480mL half filled

**Description:** 10% Neutral Buffered Formalin, 480ml prefills, 24/case

**Product No.:** PFNBF-480

**Lot No.:** 26056

**Date of Manufacture:** 10/25/2021

**Expiration Date:** 2023-10-01

**Materials:**

	<u>Part No.</u>	<u>Lot No.</u>	<u>Exp. Date</u>
<b>Chemical:</b>	NBF-275G	1J006	2023-10-01
		1J008	2023-10-01

**QC Testing:**

#### LOT NUMBER: 1J006

PROPERTY	SPECIFICATION	RESULTS	PASS / FAIL
<b>pH</b>	6.90-7.10 @ 25°C	7.06 @ 25°C	<b>PASS</b>
<b>Formaldehyde Concentration</b>	9.5-10.5% (v/v)	10.31% (v/v)	<b>PASS</b>
<b>Appearance</b>	Clear with No Particulate	Clear	<b>PASS</b>
<b>Labeling</b>	Correct Lot Number and Exp. Date	Correct	<b>PASS</b>

#### LOT NUMBER: 1J008

PROPERTY	SPECIFICATION	RESULTS	PASS / FAIL
<b>pH</b>	6.90-7.10 @ 25°C	7.01 @ 25°C	<b>PASS</b>
<b>Formaldehyde Concentration</b>	9.5-10.5% (v/v)	10.00% (v/v)	<b>PASS</b>
<b>Appearance</b>	Clear with No Particulate	Clear	<b>PASS</b>
<b>Labeling</b>	Correct Lot Number and Exp. Date	Correct	<b>PASS</b>

**Country of Origin:** United States of America

**Latex Free:** This product contains no latex

Created By:  
 James Osborn  
 QA Associate