

October 27, 2021

CERTIFICATE OF ANALYSIS

Product Name: 10% Neutral Buffered Formalin, 5 Gallon cube

Description: 10% Neutral Buffered Formalin, Tissue Preservative, 5 Gallon Cube with Spigot

Product No.: CUNBF-5-G

Lot No.: 26623

Date of Manufacture: 2021-10-26

Expiration Date: 2023-10-01

Materials:

	Part No.	Lot No.	Exp. Date
Chemical:	NBF-275G	1J009	2023-10-01
		1J010	2023-10-01
		1J011	2023-10-01

QC Testing:

LOT NUMBER: 1J009

PROPERTY	SPECIFICATION	RESULTS	PASS / FAIL
pH	6.90-7.10 @ 25°C	7.09 @ 25°C	PASS
Formaldehyde Concentration	9.5-10.5% (v/v)	10.00% (v/v)	PASS
Appearance	Clear with No Particulate	Clear	PASS
Labeling	Correct Lot Number and Exp. Date	Correct	PASS

LOT NUMBER: 1J010

PROPERTY	SPECIFICATION	RESULTS	PASS / FAIL
pH	6.90-7.10 @ 25°C	7.03 @ 25°C	PASS
Formaldehyde Concentration	9.5-10.5% (v/v)	9.70% (v/v)	PASS
Appearance	Clear with No Particulate	Clear	PASS
Labeling	Correct Lot Number and Exp. Date	Correct	PASS

LOT NUMBER: 1J011

PROPERTY	SPECIFICATION	RESULTS	PASS / FAIL
pH	6.90-7.10 @ 25°C	7.00 @ 25°C	PASS
Formaldehyde Concentration	9.5-10.5% (v/v)	10.41% (v/v)	PASS
Appearance	Clear with No Particulate	Clear	PASS
Labeling	Correct Lot Number and Exp. Date	Correct	PASS

Country of Origin: United States of America

Latex Free: This product contains no latex

Created By:
 James Osborn
 Quality Associate