



Thursday, October 28, 2021

CERTIFICATE OF ANALYSIS

Product Name: 60mL Pre-Fill Formalin

Description: Prefilled 10% Neutral Buffered Formalin, 60mL, half-filled, Formalock closure 27/pack, 2 packs/case

Product No.: PFNBF-60

Lot No.: 26733

Date of Manufacture: 2021-10-26

Expiration Date: 2023-10-01

Materials:

	<u>Part No.</u>	<u>Lot No.</u>	<u>Exp. Date</u>
Chemical:	NBF-275F	1J003	2023-10-01
		1J008	2023-10-01

QC Testing:

Formalin:

LOT NUMBER: 1J003			
PROPERTY	SPECIFICATION	RESULTS	PASS / FAIL
pH	6.90-7.10 @ 25°C	7.04 @ 25°C	PASS
Formaldehyde Concentration	9.5-10.5% (v/v)	10.31% (v/v)	PASS
Appearance	Clear with No Particulate	Clear	PASS
Labeling	Correct Lot Number and Exp. Date	Correct	PASS
LOT NUMBER: 1J008			
PROPERTY	SPECIFICATION	RESULTS	PASS / FAIL
pH	6.90-7.10 @ 25°C	7.01 @ 25°C	PASS
Formaldehyde Concentration	9.5-10.5% (v/v)	10.00% (v/v)	PASS
Appearance	Clear with No Particulate	Clear	PASS
Labeling	Correct Lot Number and Exp. Date	Correct	PASS

Country of Origin: United States of America

Latex Free: This product contains no latex

Created By:
 Fernando Fana
 Quality Associate