

Thursday, October 28,

2021

**CERTIFICATE OF ANALYSIS**

**Product Name:** 10% Neutral Buffered Formalin, 1.25 Gallon  
**Description:** Prefilled Formalin Containers, 1.25 gallon, 1/case  
**Product No.:** PFNBF-1.25G  
**Lot No.:** 26687  
**Date of Manufacture:** 10/27/2021  
**Expiration Date:** 2023-10-01

**Materials:**

	<u>Part No.</u>	<u>Lot No.</u>	<u>Exp. Date</u>
<b>Chemical:</b>	NBF-275F	1J008	2023-10-01

**QC Testing:**

<b>LOT NUMBER: 1J008</b>			
<b>PROPERTY</b>	<b>SPECIFICATION</b>	<b>RESULTS</b>	<b>PASS / FAIL</b>
<b>pH</b>	6.90-7.10 @ 25°C	7.01 @ 25°C	<b>PASS</b>
<b>Formaldehyde Concentration</b>	9.5-10.5% (v/v)	10.00% (v/v)	<b>PASS</b>
<b>Appearance</b>	Clear with No Particulate	Clear	<b>PASS</b>
<b>Labeling</b>	Correct Lot Number and Exp. Date	Correct	<b>PASS</b>

**Country of Origin:** United States of America

**Latex Free:** This product contains no latex

Created By:  
James Osborn  
QA Associate