

11/01/2021

CERTIFICATE OF ANALYSIS

Product Name: 10% Neutral Buffered Formalin, 1.25 Gallon
Description: Prefilled Formalin Containers, 1.25 gallon, 1/case
Product No.: PFNBF-1.25G
Lot No.: 26745
Date of Manufacture: 10/29/2021
Expiration Date: 2023-10-01

Materials:

	<u>Part No.</u>	<u>Lot No.</u>	<u>Exp. Date</u>
Chemical:	NBF-275G	1J010	2023-10-01

QC Testing:

LOT NUMBER: 1J010			
PROPERTY	SPECIFICATION	RESULTS	PASS / FAIL
pH	6.90-7.10 @ 25°C	7.03 @ 25°C	PASS
Formaldehyde Concentration	9.5-10.5% (v/v)	9.70% (v/v)	PASS
Appearance	Clear with No Particulate	Clear	PASS
Labeling	Correct Lot Number and Exp. Date	Correct	PASS

Country of Origin: United States of America

Latex Free: This product contains no latex

Created By:
James Osborn
QA Associate