



701 Hemlock Road
Morgantown PA 19543
Phone: 610.524.5810 Fax: 610.901.3016
info@AzerSci.com www.AzerSci.com

Certificate of Analysis

11/30/2021

Product No.: PFNBF-1000

Description: Prefilled 10% Neutral Buffered Formalin, 1000 mL, half-filled, 12/case

Lot No.: 26327

Date of Manufacturing: 11/24/2021

Expiration Date: 11/1/2023

10% Formalin Lot #: 1K009

| PROPERTY | TEST SPECIFICATION | RESULTS |
|----------------------------|---------------------------------------------|-------------|
| PH | 6.90-7.10 @ 25°C | 7.02@ 25°C |
| Formaldehyde Concentration | 9.5-10.5% (v/v) | 10.46%(v/v) |
| Appearance | Clear, colorless, no observable particulate | Conforms |

10% Formalin Lot #: 1K010

| PROPERTY | TEST SPECIFICATION | RESULTS |
|----------------------------|---------------------------------------------|-------------|
| PH | 6.90-7.10 @ 25°C | 7.07@ 25°C |
| Formaldehyde Concentration | 9.5-10.5% (v/v) | 10.46%(v/v) |
| Appearance | Clear, colorless, no observable particulate | Conforms |

Country of Origin: United States of America

Latex Free: This product contains no latex

Created By:

James E Osborn
Quality Associate