

**CERTIFICATE OF ANALYSIS**

**Product Name:** 10% Neutral Buffered Formalin, 30 Gallon Drum  
**Description:** 10% Neutral Buffered Formalin, Tissue Preservative, 30 Gallon Plastic Drum  
**Product No.:** NBF30G  
**Lot No.:** 26503  
**Date of Manufacture:** 2021-11-03  
**Expiration Date:** 2023-11-01

**Materials:**

**Chemical:** NBF-275G                      1K002    2023-11-01

**QC Testing:**

**LOT NUMBER: 1K002**

PROPERTY	SPECIFICATION	RESULTS	PASS / FAIL
pH	6.90-7.10 @ 25°C	7.02 @ 25°C	PASS
Formaldehyde Concentration	9.5-10.5% (v/v)	10.31% (v/v)	PASS
Appearance	Clear with No Particulate	Clear	PASS
Labeling	Correct Lot Number and Exp. Date	Correct	PASS

**Country of Origin:** United States of America

**Latex Free:** This product contains no latex

Created By:  
James Osborn  
Quality Associate