

## Certificate of Analysis

**Product Name:** 10% Neutral Buffered Formalin, 55 Gallon Drum

**Description:** 10% Neutral Buffered Formalin, Tissue Preservative, 55 Gallon Plastic Drum

**Product No.:** NBF55G

**Lot No.:** 26885

**Date of Manufacture:** 2021-11-04

**Expiration Date:** 2023-11-01

**Materials:**

	<u>Part No.</u>	<u>Lot No.</u>	<u>Exp. Date</u>
<b>Chemical:</b>	NBF-275G	1K003	2023-11-01
		1K004	2023-11-01

**QC Testing:**

<b>LOT NUMBER: 1K003</b>			
PROPERTY	SPECIFICATION	RESULTS	PASS / FAIL
pH	6.90-7.10 @ 25°C	7.05 @ 25°C	PASS
Formaldehyde Concentration	9.5-10.5% (v/v)	10.31% (v/v)	PASS
Appearance	Clear with No Particulate	Clear	PASS
Labeling	Correct Lot Number and Exp. Date	Correct	PASS
<b>LOT NUMBER: 1K004</b>			
PROPERTY	SPECIFICATION	RESULTS	PASS / FAIL
pH	6.90-7.10 @ 25°C	7.08 @ 25°C	PASS
Formaldehyde Concentration	9.5-10.5% (v/v)	10.00% (v/v)	PASS
Appearance	Clear with No Particulate	Clear	PASS
Labeling	Correct Lot Number and Exp. Date	Correct	PASS

**Country of Origin:** United States of America

**Latex Free:** This product contains no latex

Created By:  
 James Osborn  
 Quality Associate