

Certificate of Analysis

Product Name: 10% Neutral Buffered Formalin, 5 Gallon cube

Description: 10% Neutral Buffered Formalin, Tissue Preservative, 5 Gallon Cube with Spigot

Product No.: CUNBF-5-G

Lot No.: 26895

Date of Manufacture: 2021-11-05

Expiration Date: 2023-11-01

Materials:

	<u>Part No.</u>	<u>Lot No.</u>	<u>Exp. Date</u>
Chemical:	NBF-275G	1K001	2023-11-01
		1K003	2023-11-01

QC Testing:

LOT NUMBER: 1K001			
PROPERTY	SPECIFICATION	RESULTS	PASS / FAIL
pH	6.90-7.10 @ 25°C	7.05 @ 25°C	PASS
Formaldehyde Concentration	9.5-10.5% (v/v)	10.46% (v/v)	PASS
Appearance	Clear with No Particulate	Clear	PASS
Labeling	Correct Lot Number and Exp. Date	Correct	PASS
LOT NUMBER: 1K003			
PROPERTY	SPECIFICATION	RESULTS	PASS / FAIL
pH	6.90-7.10 @ 25°C	7.05 @ 25°C	PASS
Formaldehyde Concentration	9.5-10.5% (v/v)	10.31% (v/v)	PASS
Appearance	Clear with No Particulate	Clear	PASS
Labeling	Correct Lot Number and Exp. Date	Correct	PASS

Country of Origin: United States of America

Latex Free: This product contains no latex

Created By:
James Osborn
QA Associate