

Certificate of Analysis

Product Name: 10% Neutral Buffered Formalin, 2.5 Gallon cube

Description: 10% Neutral Buffered Formalin, Tissue Preservative, 2.5 Gallon Cube with Spigot

Product No.: CUNBF-2.5G

Lot No.: 26797

Date of Manufacture: 2021-11-10

Expiration Date: 2023-11-01

Materials:

	<u>Part No.</u>	<u>Lot No.</u>	<u>Exp. Date</u>
Chemical:	RMC-FORMALIN	1K004	2023-11-01

QC Testing:

LOT NUMBER: 1K004			
PROPERTY	SPECIFICATION	RESULTS	PASS / FAIL
pH	6.90-7.10 @ 25°C	7.08 @ 25°C	PASS
Formaldehyde Concentration	9.5-10.5% (v/v)	10.00% (v/v)	PASS
Appearance	Clear with No Particulate	Clear	PASS
Labeling	Correct Lot Number and Exp. Date	Correct	PASS

Country of Origin: United States of America

Latex Free: This product contains no latex

Created By:
James Osborn
QA Associate