

Certificate of Analysis

Product No.: 20CUNBF5G

Description: 20% Neutral Buffered Formalin, 5 Gallon Cube with Spigot

Lot No.: 26968

Date of Manufacture: 2021-11-11

Expiration Date: 2023-11-11

	<u>Part No.</u>	<u>Lot No.</u>	<u>Exp. Date</u>
Chemical:	RMC-20NBF-275	1K201	2023-11-01

QC Testing:

20% FORMALIN LOT: 1K201		
PROPERTY	SPECIFICATION	RESULTS
pH	6.90-7.10 @ 25°C	6.97 @ 25°C
Formaldehyde Concentration	19.5-20.5% (v/v)	20.32% (v/v)
Appearance	Clear with No Particulate	CONFORMS

Country of Origin: United States of America

Latex Free: This product contains no latex

Created By:
James Osborn
Quality Associate