

Certificate of Analysis

Description: Isotonic Saline Solution, 0.85% w/v, Buffered

Product No.: ES124310L

Lot. No.: 26949

Date of Manufacture: 2021-11-11

Expiration Date: 2023-05-10

Materials:

	<u>Part No.</u>	<u>Lot No.</u>	<u>Exp. Date</u>
Chemical:	ES1243-1000	1K26948	2023-05-10

QC Testing:

LOT NUMBER: 1K26948			
PROPERTY	SPECIFICATION	RESULTS	PASS / FAIL
Appearance	Clear, colorless, no observable particulate	CONFORMS	PASS
pH	7.0-7.2	7.06	PASS
Assay	0.85% to 0.90%	0.88%	PASS
TAMC	< 10 CFU/ml	<1 CFU/ml	PASS
TYMC	< 10 CFU/ml	<1 CFU/ml	PASS

Country of Origin: United States of America

Latex Free: This product contains no latex

Created By:
James Osborn
Quality Associate