

## Certificate of Analysis

**Product No.:** F-PFZF-060

**Description:** Prefilled Zinc Formalin, custom zinc formalin, 60ml, 54/cs

**Lot No.:** 26771

**Date of Manufacture:** 2021-11-12

**Expiration Date:** 2022-08-30

	<u>Part No.</u>	<u>Lot No.</u>	<u>Exp. Date</u>
<b>Chemical:</b>	<b>22-110-683</b>	457899T 109529	2022-08-30 2023-08-12

### QC Testing:

<b>22-110-683 LOT: 457899T</b>		
PROPERTY	SPECIFICATION	RESULTS
pH	6.00-8.00 @ 25°C	6.25 @ 25°C
Formaldehyde Concentration	9.0-11.4% (v/v)	9.6% (v/v)
Zinc, ppm	1700-2700ppm	2181ppm
Appearance	Clear with No Particulate	CONFORMS
<b>22-110-683 LOT: 109529</b>		
PROPERTY	SPECIFICATION	RESULTS
pH	6.00-8.00 @ 25°C	6.26 @ 25°C
Formaldehyde Concentration	9.0-11.4% (v/v)	9.6% (v/v)
Zinc, ppm	1700-2700ppm	2175ppm
Appearance	Clear with No Particulate	CONFORMS

**Country of Origin:** United States of America

**Latex Free:** This product contains no latex

Created By:  
James Osborn  
Quality Associate