



701 Hemlock Road
Morgantown PA 19543
Phone: 610.524.5810 Fax: 610.901.3016
info@AzerSci.com www.AzerSci.com

Certificate of Analysis

12/7/2021

Product No.: B-PFNBF-20

Description: Prefilled 10% Neutral Buffered Formalin, 20mL, half-filled, Formalock closure, bulk packaging 32/pack, 24 packs/case

Lot No.: 26932

Date of Manufacturing: 12/2/2021

Expiration Date: 11/1/2023

10% Formalin Lot #: 1K011

PROPERTY	TEST SPECIFICATION	RESULTS
PH	6.90-7.10 @ 25°C	7.09@ 25°C
Formaldehyde Concentration	9.5-10.5% (v/v)	10.31%(v/v)
Appearance	Clear, colorless, no observable particulate	Conforms

10% Formalin Lot #: 1K013

PROPERTY	TEST SPECIFICATION	RESULTS
PH	6.90-7.10 @ 25°C	6.94@ 25°C
Formaldehyde Concentration	9.5-10.5% (v/v)	10.16%(v/v)
Appearance	Clear, colorless, no observable particulate	Conforms

Country of Origin: United States of America

Latex Free: This product contains no latex

Created By:

James E Osborn
Quality Associate