



701 Hemlock Road  
Morgantown PA 19543  
Phone: 610.524.5810 Fax: 610.901.3016  
info@AzerSci.com www.AzerSci.com

# Certificate of Analysis

12/14/2021

**Product No.:** B-PFNBF120

**Description:** Prefilled 10% Neutral Buffered Formalin, 120mL container, half-filled, Formalock closure, bulk packaging, 24/pack, 8 packs/case

**Lot No.:** 26994

**Date of Manufacturing:** 12/13/2021

**Expiration Date:** 12/1/2023

## 10% Formalin Lot #: 1L002

PROPERTY	TEST SPECIFICATION	RESULTS
PH	6.90-7.10 @ 25°C	6.98@ 25°C
Formaldehyde Concentration	9.5-10.5% (v/v)	10.00%(v/v)
Appearance	Clear, colorless, no observable particulate	Conforms

**Country of Origin:** United States of America

**Latex Free:** This product contains no latex

**Created By:**

James E Osborn  
Quality Associate