



# Certificate of Analysis

12/29/2021

701 Hemlock Road  
Morgantown PA 19543  
Phone: 610.524.5810 Fax: 610.901.3016  
info@AzerSci.com www.AzerSci.com

**Product No.:** NBF-1G

**Description:** Formalin, 10% Neutral Buffered, 1 gallon, 4/case

**Lot No.:** 27290

**Date of Manufacturing:** 12/22/2021

**Expiration Date:** 12/1/2023

## 10% Formalin Lot #: 1L008

| PROPERTY                   | TEST SPECIFICATION                          | RESULTS     |
|----------------------------|---|-------------|
| PH                         | 6.90-7.10 @ 25°C                            | 6.95 @ 25°C |
| Formaldehyde Concentration | 9.5-10.5% (v/v)                             | 10.31%(v/v) |
| Appearance                 | Clear, colorless, no observable particulate | Conforms    |

**Country of Origin:** United States of America

**Latex Free:** This product contains no latex

**Created By:**

Fernando Fana  
Quality Associate