



701 Hemlock Road
Morgantown PA 19543
Phone: 610.524.5810 Fax: 610.901.3016
info@AzerSci.com www.AzerSci.com

Certificate of Analysis

12/22/2021

Product No.: PFNBF-240

Description: Prefilled 10% Neutral Buffered Formalin, 240mL, half-filled 24/case

Lot No.: 27296

Date of Manufacturing: 12/21/2021

Expiration Date: 12/1/2023

10% Formalin Lot #: 1L005

PROPERTY	TEST SPECIFICATION	RESULTS
PH	6.90-7.10 @ 25°C	7.02@ 25°C
Formaldehyde Concentration	9.5-10.5% (v/v)	10.31%(v/v)
Appearance	Clear, colorless, no observable particulate	Conforms

10% Formalin Lot #: 1L006

PROPERTY	TEST SPECIFICATION	RESULTS
PH	6.90-7.10 @ 25°C	7.04@ 25°C
Formaldehyde Concentration	9.5-10.5% (v/v)	10.16%(v/v)
Appearance	Clear, colorless, no observable particulate	Conforms

Country of Origin: United States of America

Latex Free: This product contains no latex

Created By:

James E Osborn
Quality Associate