

Certificate of Analysis

Product No.: PFG-40

Description: Prefilled Formagel, 10% Buffered Formalin Gel Fixative, 40mL container half-filled, 24/pack, 2 packs/case

Lot No.: 27180

Date of Manufacture: 2021-12-10

Expiration Date: 2023-05-12

	<u>Part No.</u>	<u>Lot No.</u>	<u>Exp. Date</u>
Chemical:	ESNBF30G-GEL	1132	2023-05-12

QC Testing:

FORMAGEL LOT: 1132		
PROPERTY	SPECIFICATION	RESULTS
pH	6.90-7.10 @ 25°C	7.0@ 25°C
Formaldehyde Concentration	9.9-10.4% (v/v)	10.00% (v/v)
Appearance	Clear with No Particulate	CONFORMS

Country of Origin: United States of America

Latex Free: This product contains no latex

Created By:
James Osborn
QA Associate