

Certificate of Analysis

Product No.: F-PFZF-060

Description: Prefilled Zinc Formalin, custom zinc formalin, 60ml, 54/cs

Lot No.: 27289

Date of Manufacture: 2021-12-20

Expiration Date: 2023-08-12

	<u>Part No.</u>	<u>Lot No.</u>	<u>Exp. Date</u>
Chemical:	22-110-683	109529	2023-08-12
		111631	2023-10-01

QC Testing:

22-110-683 LOT: 109529		
PROPERTY	SPECIFICATION	RESULTS
pH	6.00-8.00 @ 25°C	6.26 @ 25°C
Formaldehyde Concentration	9.0-11.4% (v/v)	9.6% (v/v)
Zinc, ppm	1700-2700ppm	2175 ppm
Appearance	Clear with No Particulate	CONFORMS
22-110-683 LOT: 111631		
PROPERTY	SPECIFICATION	RESULTS
pH	6.00-8.00 @ 25°C	6.38 @ 25°C
Formaldehyde Concentration	9.0-11.4% (v/v)	9.7% (v/v)
Zinc, ppm	1700-2700ppm	2140 ppm
Appearance	Clear with No Particulate	CONFORMS

Country of Origin: United States of America

Latex Free: This product contains no latex

Created By:
James Osborn
Quality Associate