

Certificate of Analysis

Product No.: PFZF-060

Description: Prefilled Zinc Formalin, 60mL, half-filled 27/pack, 2 packs/case

Lot No.: 27070

Date of Manufacture: 2021-12-28

Expiration Date: 2023-08-06

	<u>Part No.</u>	<u>Lot No.</u>	<u>Exp. Date</u>
Chemical:	RMC-CUNZF-5-G	1218	2023-08-06
		1321	2023-11-17

QC Testing:

LOT: 1218		
PROPERTY	SPECIFICATION	RESULTS
pH	5.6-6.2 @ 25°C	6.2 @ 25°C
Formaldehyde Concentration	9.9-10.5% (v/v)	10.2% (v/v)
Appearance	Clear, Colorless Solution	CONFORMS
LOT: 1321		
PROPERTY	SPECIFICATION	RESULTS
pH	5.6-6.2 @ 25°C	6.2 @ 25°C
Formaldehyde Concentration	9.9-10.5% (v/v)	10.5% (v/v)
Appearance	Clear, Colorless Solution	CONFORMS

Country of Origin: United States of America

Latex Free: This product contains no latex

Created By:
 James Osborn
 Quality Associate