

Certificate of Analysis

Product No.: ES772

Description: Formalin Concentrate, buffered 1 gallon, 4/case

Lot No.: 27806

Date of Manufacture: 2022-02-04

Expiration Date: 2023-04-01

	<u>Part No.</u>	<u>Lot No.</u>	<u>Exp. Date</u>
Chemical:	RMC-NBFCON	1D401	2023-04-01

QC Testing:

RMC-NBFCON LOT: 1D401			
PROPERTY	SPECIFICATION	RESULTS	PASS
pH	6.90-7.10 @ 25°C	6.96 @ 25°C	✓
Formaldehyde Concentration (diluted)	9.50-10.50% (v/v)	10.46% (v/v)	✓
Formaldehyde Concentration (as manufactured)	47.50-52.50% (v/v)	52.30% (v/v)	✓
Appearance	Clear with No Particulates	Clear	✓
Labeling	Correct Lot Number and Exp. Date	Correct	✓

Country of Origin: United States of America

Latex Free: This product contains no latex

Created By:
James Osborn
QA Associate