

Certificate of Analysis

Product Name: Azer Sterile Saline

Description: Prefilled saline, 3 ml fill in 10 ml vial, sterile, 100/case

Product No.: ES45045

Lot. No.: 27859

Date of Manufacture: 2022-02-07

Expiration Date: 2023-01-24

Product Specifications:

Ingredients: Sodium Chloride, USP Grade
Deionized Water, NCCLS Type I Reagent Grade

Packaging: Product undergoes a two-phase sub-micron filtration process as well as UV Sterilization before filling process begins.
Packaged in 10 mL vials (3 mL fill)

Sterilization: Electron Beam Sterilization (19.50 - 32.30 kGy Surface Dose)

Materials:

| | <u>Part No.</u> | <u>Lot No.</u> | <u>Exp. Date</u> |
|-----------------------|-----------------|----------------|------------------|
| Container/Cap: | 4004-1547 | 200849 | N/A |
| Chemical: | ES124420L | 27738 | 2023-07-24 |

QC Testing:

| LOT NUMBER: 27859 | | | |
|--------------------------|---|----------------|--------------------|
| ATTRIBUTE | SPECIFICATION | RESULTS | PASS / FAIL |
| Appearance | Clear, colorless, no observable particulate | Pass | PASS |
| pH | 4.5 – 7.0 | 5.09 | PASS |
| Assay | 0.85% to 0.90% | 0.89% | PASS |
| Sterilization | 19.50 - 32.30 kGy Surface Dose | 27.30 | PASS |

Country of Origin: United States of America

Latex Free: This product contains no latex

Created By:

James Osborn
Quality Assurance Associate