

Certificate of Analysis

Product No.: F-PFZF-060

Description: Prefilled Zinc Formalin, custom zinc formalin, 60ml, 54/cs

Lot No.: 28269

Date of Manufacture: 2022-03-08

Expiration Date: 2023-10-01

| | <u>Part No.</u> | <u>Lot No.</u> | <u>Exp. Date</u> |
|-----------|-----------------|----------------|------------------|
| Chemical: | 22-110-683 | 111631 | 2023-10-01 |
| | | 114792 | 2024-02-04 |

QC Testing:

| 22-110-683 LOT: 111631 | | |
|-------------------------------|---------------------------|-------------|
| PROPERTY | SPECIFICATION | RESULTS |
| pH | 6.00-8.00 @ 25°C | 6.38 @ 25°C |
| Formaldehyde Concentration | 9.0-11.4% (v/v) | 9.7% (v/v) |
| Zinc, ppm | 1700-2700ppm | 2140 ppm |
| Appearance | Clear with No Particulate | CONFORMS |
| 22-110-683 LOT: 114792 | | |
| PROPERTY | SPECIFICATION | RESULTS |
| pH | 6.00-8.00 @ 25°C | 6.37 @ 25°C |
| Formaldehyde Concentration | 9.0-11.4% (v/v) | 9.5% (v/v) |
| Zinc, ppm | 1700-2700ppm | 2325 ppm |
| Appearance | Clear with No Particulate | CONFORMS |

Country of Origin: United States of America

Latex Free: This product contains no latex

Created By:
James Osborn
Quality Associate