

## Certificate of Analysis

**Product No.:** PFZF-090

**Description:** Prefilled Zinc Formalin, 90mL, half-filled 12/pack, 4 packs/case

**Lot No.:** 27786

**Date of Manufacture:** 2022-03-24

**Expiration Date:** 2023-11-17

	<u>Part No.</u>	<u>Lot No.</u>	<u>Exp. Date</u>
<b>Chemical:</b>	RMC-CUNZF-5-G	1321	2023-11-17

### QC Testing:

<b>LOT: 1321</b>		
<b>PROPERTY</b>	<b>SPECIFICATION</b>	<b>RESULTS</b>
<b>pH</b>	5.6-6.2 @ 25°C	6.2 @ 25°C
<b>Formaldehyde Concentration</b>	9.9-10.5% (v/v)	10.5% (v/v)
<b>Appearance</b>	Clear, Colorless Solution	CONFORMS

**Country of Origin:** United States of America

**Latex Free:** This product contains no latex

Created By:  
James Osborn  
Quality Associate