

03/28/2022

**CERTIFICATE OF ANALYSIS**

**Product Name:** Zinc Formalin, 120mL, half-filled

**Description:** Prefilled Zinc Formalin, 120mL, half-filled 24/pack, 2 packs/case

**Product No.:** PFZF-120

**Lot No.:** 27781

**Date of Manufacture:** 2022-03-25

**Expiration Date:** 2023-11-17

**Materials:**

	<u>Part No.</u>	<u>Lot No.</u>	<u>Exp. Date</u>
<b>Chemical:</b>	CUNZF-5-G	1321	2023-11-17

**QC Testing:**

<b>LOT NUMBER: 1321</b>			
<b>PROPERTY</b>	<b>SPECIFICATION</b>	<b>RESULTS</b>	<b>PASS / FAIL</b>
<b>pH</b>	5.6-6.2 @ 25°C	6.2 @ 25°C	<b>PASS</b>
<b>Formaldehyde Concentration</b>	9.9-10.5% (v/v)	10.5% (v/v)	<b>PASS</b>
<b>Appearance</b>	Clear, colorless solution	Clear	<b>PASS</b>
<b>Physical Inspection</b>	Checked for correct label copy, lot number, expiration date, fill volume, absence of particulate matter and overall aesthetics.	Pass	<b>PASS</b>

**Country of Origin:** United States of America

**Latex Free:** This product contains no latex

Created By:  
 James Osborn  
 Quality Associate