



701 Hemlock Road
Morgantown PA 19543
Phone: 610.524.5810 Fax: 610.901.3016
info@AzerSci.com www.AzerSci.com

Certificate of Analysis

1/31/2022

Product No.: PFNBF-120TE

Description: Prefilled 10% Neutral Buffered Formalin, with tamper evident label, 120mL, half-filled, Formalock closure 24/pack, 2 packs/case

Lot No.: 27412

Date of Manufacturing: 1/28/2022

Expiration Date: 1/1/2024

10% Formalin Lot #: 2A006

PROPERTY	TEST SPECIFICATION	RESULTS
PH	6.90-7.10 @ 25°C	7 @ 25°C
Formaldehyde Concentration	9.5-10.5% (v/v)	10.00%(v/v)
Appearance	Clear, colorless, no observable particulate	Conforms

10% Formalin Lot #: 2A008

PROPERTY	TEST SPECIFICATION	RESULTS
PH	6.90-7.10 @ 25°C	6.97 @ 25°C
Formaldehyde Concentration	9.5-10.5% (v/v)	10.46%(v/v)
Appearance	Clear, colorless, no observable particulate	Conforms

Country of Origin: United States of America

Latex Free: This product contains no latex

Created By:

James E Osborn
Quality Associate