



Certificate of Analysis

1/4/2022

701 Hemlock Road
Morgantown PA 19543
Phone: 610.524.5810 Fax: 610.901.3016
info@AzerSci.com www.AzerSci.com

Product No.: B-PFNBF-20

Description: Prefilled 10% Neutral Buffered Formalin, 20mL, half-filled, Formalock closure, bulk packaging 32/pack, 24 packs/case

Lot No.: 27172

Date of Manufacturing: 12/29/2021

Expiration Date: 12/1/2023

10% Formalin Lot #: 1L009

PROPERTY	TEST SPECIFICATION	RESULTS
PH	6.90-7.10 @ 25°C	7.01 @ 25°C
Formaldehyde Concentration	9.5-10.5% (v/v)	10.16%(v/v)
Appearance	Clear, colorless, no observable particulate	Conforms

10% Formalin Lot #: 1L012

PROPERTY	TEST SPECIFICATION	RESULTS
PH	6.90-7.10 @ 25°C	7.04 @ 25°C
Formaldehyde Concentration	9.5-10.5% (v/v)	9.85%(v/v)
Appearance	Clear, colorless, no observable particulate	Conforms

Country of Origin: United States of America

Latex Free: This product contains no latex

Created By:

James E Osborn
Quality Associate