

Product No.: PFNBF-180**Description:** Prefilled 10% Neutral Buffered Formalin, 180mL, half-filled 50/case**Lot No.:** 27343**Date of Manufacturing:** 1/3/2022**Expiration Date:** 12/1/2023**10% Formalin Lot #: 1L010**

PROPERTY	TEST SPECIFICATION	RESULTS
PH	6.90-7.10 @ 25°C	7.05 @ 25°C
Formaldehyde Concentration	9.5-10.5% (v/v)	10.16%(v/v)
Appearance	Clear, colorless, no observable particulate	Conforms

10% Formalin Lot #: 1L011

PROPERTY	TEST SPECIFICATION	RESULTS
PH	6.90-7.10 @ 25°C	7.02 @ 25°C
Formaldehyde Concentration	9.5-10.5% (v/v)	10.16%(v/v)
Appearance	Clear, colorless, no observable particulate	Conforms

Country of Origin: United States of America**Latex Free:** This product contains no latex**Created By:**James E Osborn
Quality Associate