

Product No.: NBF-5G**Description:** Formalin, 10% Neutral Buffered, 5 gallon cube with spigot, Each**Lot No.:** 27346**Date of Manufacturing:** 1/12/2022**Expiration Date:** 1/1/2024**10% Formalin Lot #: 2A001**

PROPERTY	TEST SPECIFICATION	RESULTS
PH	6.90-7.10 @ 25°C	6.98@ 25°C
Formaldehyde Concentration	9.5-10.5% (v/v)	10.16%(v/v)
Appearance	Clear, colorless, no observable particulate	Conforms

10% Formalin Lot #: 2A002

PROPERTY	TEST SPECIFICATION	RESULTS
PH	6.90-7.10 @ 25°C	6.93@ 25°C
Formaldehyde Concentration	9.5-10.5% (v/v)	10.31%(v/v)
Appearance	Clear, colorless, no observable particulate	Conforms

10% Formalin Lot #: 2A003

PROPERTY	TEST SPECIFICATION	RESULTS
PH	6.90-7.10 @ 25°C	7.04@ 25°C
Formaldehyde Concentration	9.5-10.5% (v/v)	10.16%(v/v)
Appearance	Clear, colorless, no observable particulate	Conforms

Country of Origin: United States of America

Latex Free: This product contains no latex

Created By:

James E Osborn
Quality Associate